



GFRA3 Antibody, HRP conjugated

Product Code	CSB-PA009381LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60609
Immunogen	Recombinant Human GDNF family receptor alpha-3 protein (32-236AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GDNF family receptor alpha-3 (GDNF receptor alpha-3) (GDNFR-alpha-3) (GFR-alpha-3), GFRA3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GFRA3