



GNA11 Antibody, HRP conjugated

Product Code	CSB-PA009583LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29992
Immunogen	Recombinant Human Guanine nucleotide-binding protein subunit alpha-11 protein (1-359AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Guanine nucleotide-binding protein subunit alpha-11 (G alpha-11) (G-protein subunit alpha-11) (Guanine nucleotide-binding protein G(y) subunit alpha), GNA11, GA11
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GNA11