



GNA11 Antibody, HRP conjugated

| | |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-PA009583LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P29992 |
| Immunogen | Recombinant Human Guanine nucleotide-binding protein subunit alpha-11 protein (1-359AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Guanine nucleotide-binding protein subunit alpha-11 (G alpha-11) (G-protein subunit alpha-11) (Guanine nucleotide-binding protein G(y) subunit alpha), GNA11, GA11 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | GNA11 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |