



GNA11 Antibody, FITC conjugated

Product Code	CSB-PA009583LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29992
Immunogen	Recombinant Human Guanine nucleotide-binding protein subunit alpha-11 protein (1-359AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Guanine nucleotide-binding protein subunit alpha-11 (G alpha-11) (G-protein subunit alpha-11) (Guanine nucleotide-binding protein G(y) subunit alpha), GNA11, GA11
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GNA11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.