



GRIN1 Antibody, FITC conjugated

Product Code	CSB-PA009911LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05586
Immunogen	Recombinant Human Glutamate receptor ionotropic, NMDA 1 protein (274-451AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate receptor ionotropic, NMDA 1 (GluN1) (Glutamate [NMDA] receptor subunit zeta-1) (N-methyl-D-aspartate receptor subunit NR1) (NMD-R1), GRIN1, NMDAR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRIN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.