



GRIN2C Antibody, Biotin conjugated

Product Code	CSB-PA009914LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14957
Immunogen	Recombinant Human Glutamate receptor ionotropic, NMDA 2C protein (908-1007AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate receptor ionotropic, NMDA 2C (GluN2C) (Glutamate [NMDA] receptor subunit epsilon-3) (N-methyl D-aspartate receptor subtype 2C) (NMDAR2C) (NR2C), GRIN2C, NMDAR2C
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRIN2C