



GRIN2C Antibody, Biotin conjugated

Product Code	CSB-PA009914LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14957
Immunogen	Recombinant Human Glutamate receptor ionotropic, NMDA 2C protein (908-1007AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate receptor ionotropic, NMDA 2C (GluN2C) (Glutamate [NMDA] receptor subunit epsilon-3) (N-methyl D-aspartate receptor subtype 2C) (NMDAR2C) (NR2C), GRIN2C, NMDAR2C
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRIN2C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.