



HPGD Antibody, HRP conjugated

Product Code	CSB-PA010702LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15428
Immunogen	Recombinant Human 15-hydroxyprostaglandin dehydrogenase [NAD(+)] protein (14-130AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	15-hydroxyprostaglandin dehydrogenase [NAD(+)] (15-PGDH) (EC 1.1.1.141) (Prostaglandin dehydrogenase 1) (Short chain dehydrogenase/reductase family 36C member 1), HPGD, PGDH1 SDR36C1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	HPGD