



IFNW1 Antibody, HRP conjugated

Product Code	CSB-PA011061LB01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07352
Immunogen	Recombinant Bovine Interferon omega-1 protein (24-195AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon omega-1 (IFN-omega-c1) (Interferon alpha-II-1), IFNW1, IFNW
Immunogen Species	Bos taurus (Bovine)
Research Area	Others
Target Names	IFNW1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.