



Irf3 Antibody

Product Code	CSB-PA011818YA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P70671
Immunogen	Recombinant Mouse Interferon regulatory factor 3 protein (1-419AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon regulatory factor 3 (IRF-3), Irf3
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Irf3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.