



Irf3 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA011818YB01MO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P70671 |
| Immunogen | Recombinant Mouse Interferon regulatory factor 3 protein (1-419AA) |
| Raised In | Rabbit |
| Species Reactivity | Mouse |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Interferon regulatory factor 3 (IRF-3), Irf3 |
| Immunogen Species | Mus musculus (Mouse) |
| Research Area | Immunology |
| Target Names | Irf3 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |