



# IRF7 Antibody, Biotin conjugated

Product Code	CSB-PA011822LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92985
Immunogen	Recombinant Human Interferon regulatory factor 7 protein (131-289AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon regulatory factor 7 (IRF-7), IRF7
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IRF7