



KLRK1 Antibody, HRP conjugated

Product Code	CSB-PA012474LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26718
Immunogen	Recombinant Human NKG2-D type II integral membrane protein (73-216AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NKG2-D type II integral membrane protein (Killer cell lectin-like receptor subfamily K member 1) (NK cell receptor D) (NKG2-D-activating NK receptor) (CD antigen CD314), KLRK1, D12S2489E NKG2D
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	KLRK1