



KPNA4 Antibody, Biotin conjugated

Product Code	CSB-PA012486LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00629
Immunogen	Recombinant Human Importin subunit alpha-3 protein (51-239AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Importin subunit alpha-3 (Importin alpha Q1) (Qip1) (Karyopherin subunit alpha-4), KPNA4, QIP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KPNA4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.