



MBD3 Antibody, FITC conjugated

Product Code	CSB-PA013532LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95983
Immunogen	Recombinant Human Methyl-CpG-binding domain protein 3 protein (79-209AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Methyl-CpG-binding domain protein 3 (Methyl-CpG-binding protein MBD3), MBD3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MBD3