



Mut Antibody, Biotin conjugated

Product Code	CSB-PA015243YD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16332
Immunogen	Recombinant Mouse Methylmalonyl-CoA mutase, mitochondrial protein (31-748AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Methylmalonyl-CoA mutase, mitochondrial (MCM) (EC 5.4.99.2) (Methylmalonyl-CoA isomerase), Mut
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Mut
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.