





MYL9 Antibody

Product Code	CSB-PA015318LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24844
Immunogen	Recombinant Human Myosin regulatory light polypeptide 9 protein (2-172AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin regulatory light polypeptide 9 (20 kDa myosin light chain) (LC20) (MLC-2C) (Myosin RLC) (Myosin regulatory light chain 2, smooth muscle isoform) (Myosin regulatory light chain 9) (Myosin regulatory light chain MRLC1), MYL9, MLC2 MRLC1 MYRL2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYL9

Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA015318LA01HU at dilution of 1:100



Immunofluorescent analysis of MCF-7 cells using CSB-PA015318LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)