



PQBP1 Antibody, FITC conjugated

Product Code	CSB-PA018594LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60828
Immunogen	Recombinant Human Polyglutamine-binding protein 1 protein (2-265AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polyglutamine-binding protein 1 (PQBP-1) (38 kDa nuclear protein containing a WW domain) (Npw38) (Polyglutamine tract-binding protein 1), PQBP1, NPW38
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PQBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.