



RASGRP1 Antibody, HRP conjugated

Product Code	CSB-PA019362LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95267
Immunogen	Recombinant Human RAS guanyl-releasing protein 1 protein (594-716AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RAS guanyl-releasing protein 1 (Calcium and DAG-regulated guanine nucleotide exchange factor II) (Ca/DAG-GEFII) (Ras guanyl-releasing protein), RASGRP1, RASGRP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RASGRP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.