



SCD Antibody, Biotin conjugated

Product Code	CSB-PA020802LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00767
Immunogen	Recombinant Human Acyl-CoA desaturase protein (10-70AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acyl-CoA desaturase (EC 1.14.19.1) (Delta(9)-desaturase) (Delta-9 desaturase) (Fatty acid desaturase) (Stearoyl-CoA desaturase) (hSCD1), SCD
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SCD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.