





SHC1 Antibody

Product Code	CSB-PA021253LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29353
Immunogen	Recombinant Human SHC-transforming protein 1 protein (207-370AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SHC-transforming protein 1 (SHC-transforming protein 3) (SHC-transforming protein A) (Src homology 2 domain-containing-transforming protein C1) (SH2 domain protein C1), SHC1, SHC SHCA
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SHC1

Image



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA021253LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA021253LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.