



SHC1 Antibody, Biotin conjugated

Product Code	CSB-PA021253LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29353
Immunogen	Recombinant Human SHC-transforming protein 1 protein (207-370AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SHC-transforming protein 1 (SHC-transforming protein 3) (SHC-transforming protein A) (Src homology 2 domain-containing-transforming protein C1) (SH2 domain protein C1), SHC1, SHC SHCA
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SHC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.