



TECPR2 Antibody, HRP conjugated

Product Code	CSB-PA023369LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15040
Immunogen	Recombinant Human Tectonin beta-propeller repeat-containing protein 2 protein (516-783AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tectonin beta-propeller repeat-containing protein 2 (WD repeat-containing protein KIAA0329/KIAA0297), TECPR2, KIAA0297 KIAA0329
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TECPR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.