




TNFAIP3 Antibody

Product Code	CSB-PA023958LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21580
Immunogen	Recombinant Human Tumor necrosis factor alpha-induced protein 3 protein (97-183AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:2000-1:10000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor alpha-induced protein 3 (TNF alpha-induced protein 3) (EC 2.3.2.-) (EC 3.4.19.12) (OTU domain-containing protein 7C) (Putative DNA-binding protein A20) (Zinc finger protein A20) [Cleaved into: A20p50; A20p37], TNFAIP3, OTUD7C
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TNFAIP3

Image



Western Blot

Positive WB detected in: RAW264.7 whole cell lysate, HeLa whole cell lysate, A549 whole cell lysate


All lanes: TNFAIP3 antibody at 2.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 90 kDa

Observed band size: 90 kDa



Immunofluorescent analysis of HeLa cells using CSB-PA023958LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)