



TNFAIP3 Antibody, Biotin conjugated

Product Code	CSB-PA023958LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21580
Immunogen	Recombinant Human Tumor necrosis factor alpha-induced protein 3 protein (97-183AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor alpha-induced protein 3 (TNF alpha-induced protein 3) (EC 2.3.2.-) (EC 3.4.19.12) (OTU domain-containing protein 7C) (Putative DNA-binding protein A20) (Zinc finger protein A20) [Cleaved into: A20p50; A20p37], TNFAIP3, OTUD7C
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TNFAIP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.