




TP53I11 Antibody

Product Code	CSB-PA024082LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14683
Immunogen	Recombinant Human Tumor protein p53-inducible protein 11 protein (1-41AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor protein p53-inducible protein 11 (p53-induced gene 11 protein), TP53I11, PIG11
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TP53I11

Image



Immunofluorescent analysis of HeLa cells using CSB-PA024082LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)