



TP53I11 Antibody, HRP conjugated

Product Code	CSB-PA024082LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14683
Immunogen	Recombinant Human Tumor protein p53-inducible protein 11 protein (1-41AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor protein p53-inducible protein 11 (p53-induced gene 11 protein), TP53I11, PIG11
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TP53I11