



TP53I11 Antibody, Biotin conjugated

Product Code	CSB-PA024082LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14683
Immunogen	Recombinant Human Tumor protein p53-inducible protein 11 protein (1-41AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor protein p53-inducible protein 11 (p53-induced gene 11 protein), TP53I11, PIG11
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TP53I11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.