



TPH1 Antibody, Biotin conjugated

Product Code	CSB-PA024100LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17752
Immunogen	Recombinant Human Tryptophan 5-hydroxylase 1 protein (257-405AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tryptophan 5-hydroxylase 1 (EC 1.14.16.4) (Tryptophan 5-monooxygenase 1), TPH1, TPH TPRH TRPH
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TPH1