




VGF Antibody

Product Code	CSB-PA025847LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15240
Immunogen	Recombinant Human Neurosecretory protein VGF protein (23-206AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurosecretory protein VGF [Cleaved into: Neuroendocrine regulatory peptide-1 (NERP-1); Neuroendocrine regulatory peptide-2 (NERP-2); Antimicrobial peptide VGF[554-577]], VGF
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VGF

Image



Immunohistochemistry of paraffin-embedded human melanoma using CSB-PA025847LA01HU at dilution of 1:100