



VGF Antibody, Biotin conjugated

Product Code	CSB-PA025847LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15240
Immunogen	Recombinant Human Neurosecretory protein VGF protein (23-206AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurosecretory protein VGF [Cleaved into: Neuroendocrine regulatory peptide-1 (NERP-1); Neuroendocrine regulatory peptide-2 (NERP-2); Antimicrobial peptide VGF[554-577]], VGF
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VGF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.