



FEZ1 Antibody, HRP conjugated

Product Code	CSB-PA03389B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99689
Immunogen	Recombinant Human Fasciculation and elongation protein zeta-1 protein (12-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fasciculation and elongation protein zeta-1 (Zygin I) (Zygin-1), FEZ1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	FEZ1