



CDX2 Antibody

Product Code	CSB-PA11699A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99626
Immunogen	Recombinant Human Homeobox protein CDX-2 protein (12-116AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein CDX-2 (CDX-3) (Caudal-type homeobox protein 2), CDX2, CDX3
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	CDX2