



CDX2 Antibody, Biotin conjugated

Product Code	CSB-PA11699D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99626
Immunogen	Recombinant Human Homeobox protein CDX-2 protein (12-116AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein CDX-2 (CDX-3) (Caudal-type homeobox protein 2), CDX2, CDX3
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	CDX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.