



PTPN11 Antibody, HRP conjugated

Product Code	CSB-PA11779B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06124
Immunogen	Recombinant Human Tyrosine-protein phosphatase non-receptor type 11 protein (246-520AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein phosphatase non-receptor type 11 (EC 3.1.3.48) (Protein-tyrosine phosphatase 1D) (PTP-1D) (Protein-tyrosine phosphatase 2C) (PTP-2C) (SH-PTP2) (SHP-2) (Shp2) (SH-PTP3), PTPN11, PTP2C SHPTP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PTPN11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.