



WEE1 Antibody, HRP conjugated

Product Code	CSB-PA16059B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30291
Immunogen	Recombinant Human Wee1-like protein kinase protein (45-177AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Wee1-like protein kinase (WEE1hu) (EC 2.7.10.2) (Wee1A kinase), WEE1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	WEE1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.