



WEE1 Antibody, Biotin conjugated

Product Code	CSB-PA16059D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30291
Immunogen	Recombinant Human Wee1-like protein kinase protein (45-177AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Wee1-like protein kinase (WEE1hu) (EC 2.7.10.2) (Wee1A kinase), WEE1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	WEE1