



AQP2 Antibody, Biotin conjugated

Product Code	CSB-PA16069D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41181
Immunogen	Recombinant Human Aquaporin-2 protein (225-271AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aquaporin-2 (AQP-2) (ADH water channel) (Aquaporin-CD) (AQP-CD) (Collecting duct water channel protein) (WCH-CD) (Water channel protein for renal collecting duct), AQP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AQP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.