



AQP2 Antibody, Biotin conjugated

Product Code	CSB-PA16069D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41181
Immunogen	Recombinant Human Aquaporin-2 protein (225-271AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aquaporin-2 (AQP-2) (ADH water channel) (Aquaporin-CD) (AQP-CD) (Collecting duct water channel protein) (WCH-CD) (Water channel protein for renal collecting duct), AQP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AQP2