



MBP Antibody, Biotin conjugated

Product Code	CSB-PA20287D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	G7PWA0
Immunogen	Recombinant Macaca fascicularis Myelin basic protein (1-304AA)
Raised In	Rabbit
Species Reactivity	Macaca fascicularis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myelin basic protein, MBP, EGM_08965
Immunogen Species	Macaca fascicularis
Research Area	Others
Target Names	MBP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.