



IpxC Antibody, FITC conjugated

Product Code	CSB-PA342468LC01EZC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47205
Immunogen	Recombinant Pseudomonas aeruginosa UDP-3-O-acyl-N-acetylglucosamine deacetylase protein (1-303AA)
Raised In	Rabbit
Species Reactivity	Pseudomonas aeruginosa
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UDP-3-O-acyl-N-acetylglucosamine deacetylase (UDP-3-O-acyl-GlcNAc deacetylase) (EC 3.5.1.108) (UDP-3-O-[R-3-hydroxymyristoyl]-N-acetylglucosamine deacetylase), IpxC, envA
Immunogen Species	Pseudomonas aeruginosa (strain ATCC 15692 / DSM 22644 / CIP 104116 / JCM 14847 / LMG 12228 / 1C / PRS 101 / PAO1)
Research Area	Others
Target Names	IpxC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.