



# IpxC Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA342468LD01EZ X
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P47205
<b>Immunogen</b>	Recombinant Pseudomonas aeruginosa UDP-3-O-acyl-N-acetylglucosamine deacetylase protein (1-303AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Pseudomonas aeruginosa
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	UDP-3-O-acyl-N-acetylglucosamine deacetylase (UDP-3-O-acyl-GlcNAc deacetylase) (EC 3.5.1.108) (UDP-3-O-[R-3-hydroxymyristoyl]-N-acetylglucosamine deacetylase), IpxC, envA
<b>Immunogen Species</b>	Pseudomonas aeruginosa (strain ATCC 15692 / DSM 22644 / CIP 104116 / JCM 14847 / LMG 12228 / 1C / PRS 101 / PAO1)
<b>Research Area</b>	Others
<b>Target Names</b>	IpxC
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.