



VAC14 Antibody, HRP conjugated

Product Code	CSB-PA600170LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08AM6
Immunogen	Recombinant Human Protein VAC14 homolog protein (333-523AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein VAC14 homolog (Tax1-binding protein 2), VAC14, TAX1BP2 TRX
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	VAC14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.