



FKBP8 Antibody, Biotin conjugated

Product Code	CSB-PA617913LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14318
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase FKBP8 protein (58-190AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidyl-prolyl cis-trans isomerase FKBP8 (PPIase FKBP8) (EC 5.2.1.8) (38 kDa FK506-binding protein) (38 kDa FKBP) (FKBP-38) (hFKBP38) (FK506-binding protein 8) (FKBP-8) (FKBPR38) (Rotamase), FKBP8, FKBP38
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FKBP8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.