



LRP8 Antibody, FITC conjugated

Product Code	CSB-PA623900LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14114
Immunogen	Recombinant Human Low-density lipoprotein receptor-related protein 8 protein (131-332AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Low-density lipoprotein receptor-related protein 8 (LRP-8) (Apolipoprotein E receptor 2), LRP8, APOER2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	LRP8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.