



PDS5A Antibody, HRP conjugated

Product Code	CSB-PA643713LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q29RF7
Immunogen	Recombinant Human Sister chromatid cohesion protein PDS5 homolog A protein (2-134AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sister chromatid cohesion protein PDS5 homolog A (Cell proliferation-inducing gene 54 protein) (Sister chromatid cohesion protein 112) (SCC-112), PDS5A, KIAA0648 PDS5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PDS5A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.