



ARHGAP11B Antibody

Product Code	CSB-PA664668LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3KRB8
Immunogen	Recombinant Human Rho GTPase-activating protein 11B protein (99-231AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho GTPase-activating protein 11B (Protein FAM7B1) (Rho-type GTPase-activating protein 11B), ARHGAP11B, FAM7B1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARHGAP11B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.