



# LOLPIB Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA670318HC01LQK
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q40240
<b>Immunogen</b>	Recombinant Lolium perenne Major pollen allergen Lol p 5a protein (26-307AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Lolium perenne
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Major pollen allergen Lol p 5a (Allergen Lol p Ib) (Allergen Lol p Va) (allergen Lol p 5a), LOLPIB
<b>Immunogen Species</b>	Lolium perenne (Perennial ryegrass)
<b>Research Area</b>	Others
<b>Target Names</b>	LOLPIB
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.