



LOLPIB Antibody, Biotin conjugated

Product Code	CSB-PA670318HD01LQK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q40240
Immunogen	Recombinant Lolium perenne Major pollen allergen Lol p 5a protein (26-307AA)
Raised In	Rabbit
Species Reactivity	Lolium perenne
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major pollen allergen Lol p 5a (Allergen Lol p Ib) (Allergen Lol p Va) (allergen Lol p 5a), LOLPIB
Immunogen Species	Lolium perenne (Perennial ryegrass)
Research Area	Others
Target Names	LOLPIB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.