



HDT2 Antibody

Product Code	CSB-PA684969YA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q56WH4
Immunogen	Recombinant Arabidopsis thaliana Histone deacetylase HDT2 protein (1-306AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histone deacetylase HDT2 (HD-tuins protein 2) (Histone deacetylase 2b), HDT2, HD2 HD2B HDA4
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Research Area	Others
Target Names	HDT2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.